## PHOTOGEL-ULTRA SONO GEL

**Product Code: PG** 

Packing : 250gms, 5KG

pH : 7.00±0.5
 Viscosity : 25000 cps

- Water soluble contact medium for Ultrasound Transmission
  & for ECG, EEG, EMG & Similar application
- A long lasting, sparkling clear, hypoallergenic aqueous Gel
- Produce under vacuum restrcting air bubbles.
- Salt free, non corrosive, non irritant & non staining
- Use of PHOTOGEL do not have any damaging implication to the probes / tranducer
- PHOTOGEL could be of use to Radiologists, Obsetric, Gynaecologists, Cardiologists & Physiotherapists